
Auto Suture®
ENDO GIA® UNIVERSAL
 SINGLE USE STAPLER and ENDO GIA
 ROTICULATOR® 30 / 45 / 60-2.0, 2.5, 3.5 and 4.8
 SINGLE USE LOADING UNITS

Auto Suture®
ENDO GIA® UNIVERSAL
 SINGLE USE STAPLER and ENDO GIA
 ROTICULATOR® 30 / 45 / 60-2.0, 2.5, 3.5 and 4.8
 SINGLE USE LOADING UNITS

ENGLISH

**BEFORE USING PRODUCT,
 READ THE FOLLOWING INFORMATION THOROUGHLY.**

IMPORTANT!

This booklet is designed to assist in using the Auto Suture® ENDO GIA® UNIVERSAL stapler and single use loading units. It is not a reference to surgical endoscopic techniques. For information on endoscopic procedures, techniques, complications and hazards, please see the book: *Surgical Laparoscopy Update* Karl A. Zucker, Ed. Quality Medical Publishing, Inc., St. Louis MO, 1993.

CAUTION: Federal (USA) law restricts this device to sale and use by, or on the order of, a physician.

INDICATIONS

The Auto Suture® ENDO GIA® UNIVERSAL stapler has applications in abdominal, gynecologic, pediatric and thoracic surgery for resection, transection and creation of anastomosis.

EFFECTS

The Auto Suture® ENDO GIA® UNIVERSAL stapler places two, triple-staggered rows of titanium staples and simultaneously divides the tissue between the two, triple-staggered rows. The size of the staples is determined by the selection of the 2.0 mm, 2.5 mm, 3.5 mm or 4.8 mm single use loading unit (SULU). The ENDO GIA® UNIVERSAL stapler will accommodate any of the single use loading unit sizes that are available in the 30 mm, 45 mm and 60 mm lines. The ENDO GIA® UNIVERSAL stapler with either the 2.0, 2.5, 3.5 single use loading unit is designed for introduction and use through a 12 mm trocar sleeve, or larger, with the use of a converter.

When using the ENDO GIA® UNIVERSAL stapler with the 4.8 single use loading unit, it must be inserted into a 15 mm trocar sleeve.

NOTE: Each instrument can accommodate the 30-2.0, 2.5, 3.5 / 45-2.0, 2.5, 3.5, 4.8 / 60-2.5, 3.5, 4.8 mm SULU's.

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
LOADING

1. The SULU is packaged in the open position. Do not attempt to close the SULU.

CAUTION: Do not attempt to remove the shipping wedge until the SULU is loaded into the instrument.

2. Ensure that the black knobs on the instrument are pulled back completely and the articulation arm is neutral to the instrument.

A) BLACK RETURN KNOB

3. To load the ENDO GIA® UNIVERSAL stapler with the appropriate SULU, insert the pin located at the distal end of the instrument shaft into the SULU. Ensure that the white alignment loading line on the SULU aligns with the white alignment loading line on the shaft. Push the SULU in and twist clockwise 45° relative to the instrument so that the SULU will lock into

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place. At this point the lines will realign

B) WHITE ALIGNMENT LOADING LINE (SULU)

C) WHITE ALIGNMENT LOADING LINE (SHAFT)


4 Remove the shipping wedge from the SULU prior to inserting the instrument into the trocar.

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UNLOADING

1 To unload a SULU from the stapler, the articulating lever must be in the neutral position. Ensure that the jaws of the SULU are open by pulling the black return knobs back completely. Pull the UNLOAD/UNLOCK button (located on the underside of the shaft) back towards the instrument, twist the SULU counterclockwise 45° and remove the SULU from the shaft of the instrument.

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
INSTRUCTIONS FOR USE

NOTE: The jaws of the SULU must be closed prior to introducing the instrument into the trocar sleeve. To do so, squeeze the handle.

1. Insert the END O GIA® UNIVERSAL stapler into an appropriately sized trocar sleeve, or larger, with the use of a converter. The anvil must be completely visible (past the trocar sleeve) prior to opening the SULU within the body cavity.

NOTE: When using the END O GIA® UNIVERSAL instrument with a 4.8 single use loading unit, the instrument **MUST** be inserted into a 15 mm trocar. A smaller size trocar will not accept the 4.8 single use loading unit.

2. Once inside the body cavity, open the jaws of the instrument by pulling the black return knobs completely back.

DO NOT SQUEEZE THE HANDLE WHILE PULLING BACK THE BLACK RETURN KNOBS.

3. The instrument shaft rotates 360° and articulates 22° and 45° in both directions.

4. Apply the END O GIA® UNIVERSAL stapler across the tissue to be transected. Ensure that no obstructions (such as clips) are incorporated in the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples. The instrument will not cut tissue beyond the black "cut mark" indicated on the single use loading unit. More than one application of the END O GIA® UNIVERSAL stapler may be necessary for tissue exceeding the length of the SULU (30 mm, 45 mm or 50 mm).

A) CUT MARK

NOTE: Placement of tissue proximal to the tissue stops (on the SULU) may result in stapler malfunction. Any tissue extending beyond the cut mark will not be transected.

5. Close the jaws of the instrument across the tissue to be transected by squeezing the handle completely. The stapler is equipped with a safety interlock; the instrument will not fire the staples and cut tissue unless the green button is pushed.

DO NOT ATTEMPT TO OVERRIDE THE SAFETY INTERLOCK; TO DO SO WILL RENDER THE STAPLER NONOPERATIONAL.

The jaws of the instrument are atraumatic and can be repositioned on the tissue prior to firing by pulling the black return knobs completely back allowing the jaws to open.

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
6 In order to fire the instrument, push the green button. Squeeze handle sequentially until the "oval clamp cover" reaches the distal end of the cartridge slot, and the handle locks.

Sequential squeezes of the handle are required to fully fire the SULU. Number of squeezes is relative to the length of the SULU (30, 45 or 60).

FAILURE TO COMPLETELY FIRE THE SULU WILL RESULT IN AN INCOMPLETE CUT AND INCOMPLETE STAPLE FORMATION, AND MAY RESULT IN POOR HEMOSTASIS.

B) GREEN BUTTON

C) OVAL CLAMP COVER

7 Once the instrument has been completely fired, pull the black return knobs completely back, allowing the release of the tissue from the jaws. Gently remove the instrument from the tissue. The application site should be checked for hemostasis following removal of the instrument. Minor bleeding may be controlled using electrocautery or manual sutures.

8 After completely firing the SULU, close jaws of the instrument and remove the ENDO GIA[®] UNIVERSAL stapler from the body cavity to unload the single use loading unit from the instrument.

NOTE: Do not attempt to insert or remove the instrument from the trocar sleeve if the instrument is in the articulated position.

THE ENDO GIA[®] UNIVERSAL INSTRUMENT CAN BE LOADED AND FIRED UP TO 25 TIMES.

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
G) STAPLE DIAMETER

PRECAUTION

1. When using the ENDO GIA[®] UNIVERSAL instrument with a 4.8 single use loading unit, the instrument MUST be inserted into a 15 mm trocar. A smaller size trocar will not accept the 4.8 single use loading unit.

2. Always inspect the tissue thickness and select an appropriate staple size prior to application of the ENDO GIA[®] UNIVERSAL stapler.

3. Always close the jaws of the ENDO GIA[®] UNIVERSAL stapler prior to introducing and removing the stapler from the trocar sleeve.

4. After firing, always inspect the staple line for hemostasis. Minor bleeding may be controlled by electrocautery or manual sutures.

5. Placement of tissue proximal to the tissue stops (on the SULU) may result in stapler malfunction. Any tissue extending beyond the cut mark will not be transected.

6. When using the stapler more than once during a SINGLE surgical procedure, be sure to remove the empty ENDO GIA[®] UNIVERSAL single use loading unit and reload a new one. A safety interlock is provided that prevents an empty single use loading unit from being fired a second time. Do not attempt to override the safety interlock.

7. When positioning the stapler on the application site, ensure that no obstructions, such as clips, are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

8. Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. Prior to performance of any endoscopic procedures, consult the medical literature relative to techniques, complications and hazards.

9. A thorough understanding of the principles involved in laser and electrosurgical procedures is essential to avoid shock and burn hazards to both patient and operator(s), and damage to the instrument.

10. When endoscopic instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that electrical isolation or grounding is not compromised.

11. The anvil must be completely visible (past the trocar sleeve) prior to opening the SULU within the body cavity.

12. The instrument and single use loading unit are provided STERILE and are intended for use in a SINGLE procedure only. DISCARD AFTER USE. DO NOT

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
RESTERILIZE.

13. Do not attempt to load SULU while squeezing ring handle.

CONTRAINDICATIONS

1. Do not use the ENDO GIA® UNIVERSAL 20 mm staples on any tissue that compresses to less than .75 mm in thickness, on any tissue that cannot comfortably compress to 10 mm or on the aorta.
2. Do not use the ENDO GIA® UNIVERSAL 25 mm staples on any tissue that compresses to less than 10 mm in thickness, on any tissue that cannot comfortably compress to 15 mm or on the aorta.
3. Do not use the ENDO GIA® UNIVERSAL 35 mm staples on any tissue that compresses to less than 15 mm in thickness, on any tissue that cannot comfortably compress to 20 mm or on the aorta.
4. Do not use the ENDO GIA® UNIVERSAL 48 mm staples on any tissue that compresses to less than 20 mm in thickness, on any tissue that cannot comfortably compress to 25 mm or on the aorta.
5. The ENDO GIA® UNIVERSAL instrument should not be used on tissue such as liver or spleen where compressibility is such that closure of the instrument would be destructive.
6. Do not use the ENDO GIA® UNIVERSAL stapler where adequacy of hemostasis cannot be verified visually after applications.
7. These devices are provided STERILE and are intended for use in a SINGLE procedure only. DISCARD AFTER USE. DO NOT RESTERILIZE.

STORE AT ROOM TEMPERATURE. AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.

 **< 54° C** DO NOT EXPOSE
130° F TO TEMPERATURES ABOVE 130° F (54° C).

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CE 0123

*Trademark
Manufactured by: United States Surgical,
a division of Tyco Healthcare Group LP, Norwalk, CT 06856
© 2000 United States Surgical
All Rights Reserved 07/00 Made in USA Q2
European Representative: Tyco Healthcare UK Ltd, Gosport, PO13 0AS, UK

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Vertretung f r Europa: Tyco Healthcare UK Ltd, Gosport, PO13 0AS, GB
Rappresentante per l'Europa: Tyco Healthcare UK Ltd, Gosport, PO13 0AS, UK
Representante Europeo: Tyco Healthcare UK Ltd,